

Air mattress


Air mattress

**Air Mattress Cell**  
Model: AM40



Adjustable hangers to suit any kind of hospital bed frame

- Cell mattress
- 2 to 1 alternating therapy
- Adjustable hangers for most bed types
- Comfortable Nylon mattress
- High quality pump
- Weight indicator pressure adjustment setting
- 12 minutes cycling time
- Near silent low noise
- For treatment of patients at medium to high risk of pressure ulcers




**Air Mattress Bubble**  
Model: AM30



Adjustable hangers to suit any kind of hospital bed frame

- Bubble mattress
- 2 to 1 alternating therapy
- Adjustable hangers for most bed types
- Durable medical grade PVC
- Extension flags for stabilization
- High quality pump
- Weight indicator pressure adjustment setting
- 6 minutes cycling time
- Near silent low noise
- For prevention and treatment of patients at low risk of pressure ulcers



Accessories




High Quality Rubber Tube

Repairing Kit

Fuse

Accessories



High Quality Rubber Tube

Spare Cell

Fuse

Rossmax International Ltd.  
12F, No.189, Kang Chien Rd.,  
Taipei 114, Taiwan.  
Tel +886 2 2659 7888  
Fax +886 2 2659 7666  
info@rossmax.com

Rossmax Swiss GmbH  
Tramstrasse 16, 9442 Berneck  
Switzerland  
Tel +41 71 747 11 93  
Fax +41 71 747 11 91  
emea@rossmax.com

[www.rossmax.com](http://www.rossmax.com)

OBM\_OP\_Pro\_EN\_ver1545

Professional Range

PROFESSIONAL

## Sphygmomanometer

### Aneroid Sphygmomanometer Palm Type

Model: GD series



- Superior performance and quality
- Big dial face 65mm & easy to read
- Single-hand operation
- Non-stop pin reliable manometer
- Latex free cuff

GD102 - D-ring cuff with single head stethoscope



### Aneroid Sphygmomanometer

Model: GB series  
AGC



- White enamel 300mmHg non-stop pin reliable manometer
- Chrome-plate brass air release valve
- Light-weight stethoscope
- Latex free cuff

GB102 - D-ring cuff with single head stethoscope



### Digital Sphygmomanometer

Model: GA series  
Mandaus II



- Marker
- Extra large backlight LCD
- Shock and splash resistance
- Chrome-plate brass air release valve
- Light-weight stethoscope
- 360° swivel cuff clip

GA102 - D-ring cuff with single head stethoscope

• Patent No.: CN: ZL01259708.2; TW: 488518;  
GE: 20206337.2; JP: 3089650; US: 717746



## Stethoscopes

### Single Head Stethoscope

Model: EB100



- Soft eartips
- Aluminium chest piece
- Diaphragm side, detects high frequency sound
- Snap on diaphragm ring
- Latex-free tube

### Dual Head Stethoscope

Model: EB200



- Soft eartips
- Aluminium chest piece
- Non-chill bell and diaphragm ring
- Diaphragm side detects high frequency sound
- Bell side detects low frequency sound
- Snap on diaphragm ring
- Latex-free tube

### Sprague Rappaport Stethoscope

Model: EB500



5 IN 1

5 In 1 Multipurpose



Interchangeable chestpiece – adult & child diaphragm



Interchangeable chestpiece – 3 bells (1 1/4", 1", 3/4")



Proprietary valve mechanism design



Double tube configuration

### Cardiology Stethoscope

Model: EB600



Super acoustic sensitivity



Advanced stainless steel construction



Interchangeable chestpiece – adult & child diaphragm



Proprietary valve mechanism design



2 tubes in 1

## Digital Blood Pressure Monitor for Professional Use

### Digital Professional Blood Pressure Monitor


Model: AC1000f



	High performance and accuracy medical pump and valve
	Both automatic and auscultatory measurement available
	Markers for Systolic & Diastolic at auscultatory
	Different deflation rate setting
	Easily change of different size cuffs
	7 memories for automatic mode
	Battery Charging for mobility & saving energy
	Anti-splash Protective Cover
	Lifetime Calibration

300,000  
Measurements  
Guaranteed!

Hypertension is already a prevalent risk factor for development of cardiovascular disease due to increasing longevity and other contributing factors such as obesity. While the treatment of hypertension has been shown to prevent cardiovascular diseases and extend life-span, hypertension still remains inadequately treated in many places in the world. Multiple guidelines for a more effective management of hypertension have been published and strictly followed especially in clinical practice. Rossmax Professional product range is there to cater all your needs to accurately diagnose hypertension in the early stage while conforming to the highest standards and requirements.



Trolley (optional)